



Canine Rehabilitation & Fitness Center, Inc.

TEL: (619) 846-9531 FAX: (858) 430-3999 www.TsavosCanineRehab.com

PST Indications, Contraindications, Risks and Side Effects

Leased to: _____ (print Client name)

Indications

The PST VET system was developed for the treatment of musculoskeletal disorders, particularly osteoarthritis (arthrosis), including degenerative joint disease, and post-traumatic or post-operative disorders. It has also been found to have positive healing effects in OCD, non-healing wounds, tendonitis & bursitis, and meniscal tears.

PST VET uses patented signals that emulate those innate in the body, and so synergistically acts to effectively support and accelerate the body's own natural healing processes.

A course of treatment entails either 9 sessions conducted either once daily or twice daily, 30 minutes in duration, carried out on consecutive days. If conducting twice daily treatments, each treatment must be separated by at least 8 hours. If need be, a one day pause after the 6th treatment is permitted.

Treatment must not be interrupted for longer than one day (24 hours) under any circumstances!

For Portable PST Only: If you are using the portable PST, first remove the unit from the case and place the case at least 15 feet away from the unit itself (this prevents magnetic interference). To treat, simply place the PST VET applicator on the area to be treated in much the same way you would an ice or heating pad. Make sure the joint you are treating is not in a fully flexed (bent) position. Note the direction of the writing the 1st time you place it on the dog. This **MUST** be the same direction each successive treatment or you will negate (erase) the effect of some of the treatments.

PST VET can be used to treat/relieve pain in the leg, shoulder, elbow, spine, and hip, following consultation with your pet's veterinarian.

Contraindications

In general, patients with tumors and/or bacterial infections, metallic implants (pure iron), as well as pregnant dogs should not be treated. Owners of the pets with tumors and/or bacterial infections, metallic implants (pure iron), or that are pregnant should also not be present during treatment.

Risks and Side Effects

During treatment, the dog may experience a slight tickling sensation or warmth. In general, a change, or transient intensification, in pain symptoms may also occur in the initial phase of the treatment. However, do not be concerned, as this is only temporary and considered a positive prognostic response to the treatment.

Accompanying Therapy and Medication

In general, a patient's medication should not be altered during PST treatment, unless the veterinarian concerned chooses to do so. Acupuncture, chiropractic therapy, electrical stimulation or any other therapy that changes the electrical or bio-magnetic charge of the patient should be discontinued during the course of the PST treatments. If you have any question about specific therapies or treatments, please ask one of our doctors.

I have read and understand the above information regarding the PST treatment in entirety, and I acknowledge that neither my pet nor anyone administering the treatment has any of the conditions listed in the "Contraindications" statement.

(Client Signature)

(Tsavo's Canine Rehabilitation, Inc. Signature)